

Development of a population balance model for continuous twin-screw granulation in pharmaceutical manufacturing

Ashish Kumar, Krist V. Gernaey, Thomas De Beer, Ingmar Nopens

AICHE Annual Meeting 2014

Session: Agglomeration and Granulation Processes


LABORATORY OF PHARMACEUTICAL PROCESS ANALYTICAL TECHNOLOGY

FACULTY OF PHARMACEUTICAL SCIENCES

ISCOMATH, DEPARTMENT OF MATHEMATICAL MODELLING, STATISTICS AND BIOINFORMATICS

FACULTY OF BIOSCIENCE ENGINEERING

Current solid-dosage manufacturing is slow and expensive




Product collected after each unit operation

Actual processing time = days to weeks

Traditional to new granulation method


High-shear mixer



Drum



Fluidised-bed




From minutes to hours

Batch

Continuous

Twin-screw




in seconds

content for reference only, not for reuse

Easy to integrate with other unit operations of pharmaceutical manufacturing

content for reference only, not for reuse

content for reference only, not for reuse




content for reference only, not for reuse

content for reference only, not for reuse




content for reference only, not for reuse

Continuous manufacturing line



Consigma™-25 system



Both geometry and process conditions drive constitutive mechanisms



Understanding the role of screw design




content for reference only, not for reuse

Consigma™-1 system

(GEA pharma systems, Collette)

content for reference only, not for reuse

Open barrel of a twin screw granulator



content for reference only, not for reuse

content for reference only, not for reuse

content for reference only, not for reuse

content for reference only, not for reuse

content for reference only, not for reuse

Consigma™ - 1 experiments

Granulated Lactose monohydrate with distilled water

Factors:

Parameters	Low	High
Throughput	10 Kg/h	25 Kg/h
Liquid-solid ratio	4.58%	6.52%
Screw speed	500 RPM	900 RPM



1

2

3

4

5

kneading block 1



kneading block 2

Responses:

Particle characterization by Dynamic Image Analysis
(Location 1, 3, 5)

content for reference only, not for reuse


Consigma™ - 1 experiments





Throughput **High**

Liquid-solid ratio **High**


Screw speed **Low**



Consigma™ - 1 experiments




Throughput **High** Liquid-solid ratio **High** Screw speed **High**





For details see: Kumar, et al. "Experimental investigation of granule size and shape dynamics in twin-screw granulation." *I J Pharma* 475.1 (2014): 485-495.

Tracer concentration in granules produced was measured using NIR chemical imaging



API Map was used to measure RTD




For details see: Kumar, et al. "Mixing and transport during pharmaceutical twin-screw wet granulation: Experimental analysis via chemical imaging." *E J Pharma. Biopharma*. 87.2 (2014): 279-289.

content for reference only, not for reuse

Measure of the mean of the distribution

content for reference only, not for reuse



content for reference only, not for reuse

content for reference only, not for reuse

$$\tau = \frac{\int_0^{\infty} t \cdot e(t) dt}{\int_0^{\infty} e(t) dt}$$

Mean residence time, τ

content for reference only, not for reuse

content for reference only, not for reuse

Population balance models can track granule attributes

$$\frac{\partial n(t, x)}{\partial t} = \frac{Q_{in}}{\tilde{V}} n_{in}(x) - \frac{Q_{out}}{\tilde{V}} n_{out}(x)$$

GSD balance

Aggregation
term

$$+ \frac{1}{2} \int_0^x \underbrace{\beta(t, x - \varepsilon, \varepsilon)}_{\text{aggregation rate}} n(t, x - \varepsilon) n(t, \varepsilon) d\varepsilon$$

$$- n(t, x) \int_0^\infty \underbrace{\beta(t, x, \varepsilon)}_{\text{aggregation rate}} n(t, \varepsilon) d\varepsilon$$

Breakage
term

$$+ \int_0^\infty \underbrace{b(x, \varepsilon)}_{\text{breakage fun.}} \underbrace{S(\varepsilon)}_{\text{selection rate}} n(t, \varepsilon) d\varepsilon$$

$$- \underbrace{S(x)}_{\text{selection rate}} n(t, x)$$

Semi-empirical kernels

Aggregation Kernel

$$\beta(x, y) = \beta_0$$

(Constant kernel)


Breakage Kernel

$$S(y) = S_0(y)^\mu$$


$$b(x, y) = \frac{\frac{\phi\gamma x^{\gamma-1}}{y^\gamma} + \frac{(1-\phi)\alpha x^{\alpha-1}}{y^\alpha}}{\frac{\phi\gamma}{\gamma+1} + \frac{(1-\phi)\alpha}{\alpha+1}}$$

(Austin, 2002)


Cell average technique





Experimental and simulated data have a good agreement




Particle population dynamics during granulation



including effect of granulator design on granule size distribution






	Parameter	Significance	Location 1	Location 3	Location 5	End
	RMSE		16.25	2.72	2.14	4.3
Aggregation	β_0	Collision frequency (sec ⁻¹)	1052.53	1632.92	820.79	5055.40
	S_0	Selection function constant for breakage	0.07	0.99	2.40	3.31
Breakage	α	Width of fragment distributions	0.17	0.31	0.39	0.35
	γ	width of fragment distributions	22.21	28.89	2.64	200.30
	ϕ	mass content of first breakage distributions	0.02	0.00	0.02	0.01


γ , ϕ and α are dimensionless material constants.
 ϕ is the weight parameter to quantify the mass content of first breakage distributions.
 γ and α are width of the fragment distributions ϕ and $1-\phi$, respectively.
 Quadratic selection function, $S(y) = S_0 (y)^\mu$ where μ was 1/3.

content for reference only, not for reuse


Including effect of granulator design on granule size distribution



Start	1	2	3	4	5
PMSE	16.25		2.72		2.14



Including effect of granulator design on granule size distribution




1

2

3

1st kneading block





3

4

5

2nd kneading block




content for reference only, not for reuse

Throughput **High**


Liquid-solid ratio **High**

Screw speed **Low**

1 mixing zone



2 mixing zones




content for reference only, not for reuse

Throughput **High**


Liquid-solid ratio **High**

Screw speed **High**

1 mixing zone




2 mixing zones



content for reference only, not for reuse

Investigating effect of screw speed and screw configuration



High throughput, high L/S

	Low Screw Speed				High Screw Speed			
	1 mixing zone		2 mixing zones		1 mixing zone		2 mixing zones	
Zone	1-3	3-5	1-3	3-5	1-3	3-5	1-3	3-5
RMSE	2.424	2.317	2.716	3.929	1.153	3.366	8.176	3.772
R ²	0.989	0.987	0.984	0.983	0.989	0.983	0.97	0.982

content for reference only, not for reuse

Throughput **High**




Liquid-solid ratio **High**

Location 1-3



Screw speed Low

Location 3-5

mixing zones



Screw speed High



content for reference only, not for reuse

content for reference only, not for reuse




Throughput **High**

Liquid-solid ratio **High**



Location 1-3

Screw speed Low

Location 3-5



Screw speed High



mixing zones

content for reference only, not for reuse

content for reference only, not for reuse

content for reference only, not for reuse

content for reference only, not for reuse

content for reference only, not for reuse

content for reference only, not for reuse

Conclusions

Along with experimental study, an improved insight can be obtained by model-based analysis.

Wetting kinetics requires a separate explanation in the twin-screw granulation modelling.

Aggregation and breakage are most dominant phenomena in the twin-screw granulation.

Particle population dynamics and screw geometry effect can be better understood by compartmental PBM, and can ultimately be used for predictive modelling of twin-screw granulation.

content for reference only, not for reuse


Acknowledgements

content for reference only, not for reuse



DTU

content for reference only, not for reuse



content for reference only, not for reuse



content for reference only, not for reuse

content for reference only, not for reuse

Ashish.Kumar@UGent.be

content for reference only, not for reuse